Dermaroller MC

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Dermaroller MC

Medical device acc. to Regulation (EU) 2017/745 on medical devices (MDR) – MADE IN GERMANY

Please read the instructions carefully.

They contain important information on how to handle the microneedling device Dermaroller MC.

Medical purpose: Minimally invasive perforation of the dermis (microneedling)

Effect: Microneedling causes minimally invasive micro-injuries of epidermis and dermis, without causing open wounds in the classical sense. The puncturing channels close after a short time. The micro-injuries trigger the body's wound healing cascade, which, inter alia, induces an increased collagen production in the treated skin. For this reason, the method is also called "percutaneous collagen induction" (PCI) or "Collagen Induction Therapy" (CIT).

Indications:

- Medical indication:

Treatment of skin changes caused by disease and injury, especially atrophic and hypotrophic scars (e.g. acne scars) as well as mature hypertrophic on the entire body; treatment of striae distensae

- Non-medical indication:

Stimulation of collagen production with the aim of reducing the signs of ageing / wrinkles and general skin tightening on the face, neck, décolleté or body

Clinical benefit:

Improvement of scars and striae (reduction of scar depth, normalization of skin structure) as well as of wrinkles and sagging skin (reduction of wrinkle depth, normalization of skin structure).

Patient group:

In general, the Dermaroller MC is applicable to almost all skin areas of the patients and in all skin types. Exceptions are described under "Contraindications".

User and recommended needle length:

The device's use in the context of a medical indication is reserved for qualified medical professionals. Depending on the type and severity of the scar and the area of skin to be treated, the user selects the most suitable needle length up to a maximum of 2.5 mm. To minimize adverse reactions, always choose the shortest possible needle length to achieve even redness or superficial pinpoint bleeding.

For acne scars, we recommend a needle length of 1.5 mm.


In the context of a non-medical indication, non-medical professional users such as trained beauticians may also use the medical device. For this group of users, we recommend staying below a needle length of 1.5 mm and always selecting the needle length so that the treated skin area is evenly and slightly reddened, but there is no bleeding.

Precautions / Notes:

- The device is only to be used by trained professionals
- The treatment must be carried out under low-germ conditions and in hygienically suitable treatment rooms
- Sterile disposable gloves must be worn during the treatment
- Medical Device for **SINGLE USE only**. In case of reuse there would be a risk of skin infections!
- Reprocessing of the device is prohibited since this could lead to damage or quality deterioration
- Do not use the device if the sterile packaging is damaged. Dispose of the device properly.
- Store at room temperature and away from heat sources
- After the expiration date, the instrument may not be used but is to be disposed of properly.

Procedure:

- Clean the skin area to be treated, anaesthetize with a topical anaesthetic cream, then disinfect the skin prior to treatment
- Roll the same area 4-6 times in each direction (as illustrated) until the treated skin shows a uniform mild redness and pinpoint bleedings appear.
- Apply only mild pressure and do not apply lateral pressure on the device
- Don't change rolling angle while the needles are in the skin
- In skin areas with thin dermis directly over bone (e.g., periorbital) only very low pressure may be used

Post-treatment

Clean the treated skin area with sterile swabs and sterile saline solution.

Treatment cycles:

- A sustained treatment success can be expected only after repeated treatments.
- Usually 3-5 treatments are carried out, but depending on scar / wrinkle depth and structure also more treatments may be necessary.
- The interval between two treatments should be at least 4-6 weeks to make sure that the wound healing cascade is completed.

Contraindications:

- Treatment of mucous membranes or eyes (strictly contraindicated)
- Patients with or a history of keloids (strictly contraindicated)
- Treatment of skin areas with skin cancer (strictly contraindicated)
- Active herpes simplex, active skin infections, raised moles, warts, solar keratoses (except in the context of photodynamic therapy)

- Increased bleeding tendency or wound healing disorders, e.g. due to uncontrolled diabetes mellitus
- Anti-coagulant therapy, chemotherapy, radiotherapy, high doses of cortisone
- Pregnant / breastfeeding women: since there are no clinical results on this patient group, we recommend postponing the treatment
- Children / adolescents: the clinical study results were obtained predominantly in adults. Therefore, the recommendation for the treatment is "18 years and over". However, there is no evidence that treatment in adolescents and children poses an increased risk to safety.

Accessory symptoms of the treatment:

Transient redness and swelling of the treated skin, pain, microbleedings for minutes

Potential adverse effects:

Skin infections due to inadequate disinfection, inflammatory response (caused by pyrogens), transient pigmentation disorders, scarring when a tendency to keloid formation is present, granuloma after the use of active substances during the Dermaroller treatment.

Good skin hygiene before and after the procedure is advised to further minimise the risk of infection.

Prolonged sun or other UV exposure is not recommended after the procedure for at least 8 days to minimize the risk of hyperpigmentation.

Patients should use sunscreen that provides a high level of UVA and UVB protection.

The stainless steel of the Dermaroller microneedles contains ca. 8% nickel. Due to the short treatment duration, it is very unlikely that a patient develops an allergic reaction.

Information obligations:

Before the treatment, the user must inform the patient about the contraindications, side effects and possible undesirable side effects.

Both users and patients are obliged to report all *serious* incidents that have occurred in connection with the medical device to the manufacturer and the local responsible authority.

Disposal

Used Dermarollers pose a risk of injury and a potential health hazard due to contact with blood. They are classified as „sharps“ (European Waste Catalogue 18 01 01) and must be collected in separate, closed containers and then disposed of safely.

The product as well as the packaging must be disposed of in accordance with the applicable national environmental protection laws or disposal regulations.

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