NADELMODUL: NM6xx, NM6xx LS

Medizinprodukt gem. Richtlinie 93/42/EWG – MADE IN GERMANY

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Revision 2018.11.06

Needle module: NM6xx, NM6xx LS

Medical device acc. to Directive 93/42/EEC – MADE IN GERMANY

DERMAROLLER GmbH

Wilhelm-Mast-Str. 8

38304 Wolfenbuettel

Germany

Tel. +49 5331 7108480

info@dermaroller.de

Please read the instructions carefully.

They contain important information on how to use the needle module with the microneedling device eDermastamp.

Medical purpose: Minimally invasive perforation of the dermis (microneedling)

Effect: Microneedling causes minimally invasive micro-injuries of epidermis and dermis, without causing open wounds in the classical sense. The puncturing channels close after a short time. The micro-injuries trigger the body's wound healing cascade, which, inter alia, induces an increased collagen production in the treated skin. For this reason, the method is also called "percutaneous collagen induction" (PCI) or "Collagen Induction Therapy" (CIT).

Indications:

To treat and repair atrophic and hypotrophic scars including acne scars, mature hypertrophic scars resulting from injury, trauma or disease, and stretch marks (*striae distensae*).

In general, the needle module is applicable to almost all skin areas of the patients and in all skin types.


Exceptions are described under "Contraindications".

Precautions / Notes:

- The device is only to be used by trained professionals
- The treatment must be carried out under low-germ conditions and in hygienically suitable medical treatment rooms
- Sterile disposable gloves must be worn during the treatment
- Medical Device for **SINGLE USE only**. In case of reuse there would be a risk of skin infections!
- Reprocessing of the device is prohibited since this could lead to damage or quality deterioration
- Do not use if the sterile packaging is damaged
- Store at room temperature and away from heat sources
- After the expiration date, the instrument may not be used

Preparation:

- The instructions for use of the eDermastamp device must be followed!
- The hygienic precautions must be followed.
- **Mounting:** The sterile needle module is screwed on the hand piece (see Fig.):

Procedure:

- Clean the skin area to be treated, anaesthetise with a topical anaesthetic cream, then disinfect the skin prior to treatment
- Pass the needle module with light pressure in overlapping circular motions over the skin until the treated skin shows a uniform mild redness and pinpoint bleedings appear
- Do not apply lateral pressure on the device
- In skin areas with thin dermis directly over bone (eg, periorbital) only very low pressure may be used

Post-treatment:

- Clean the treated skin area with sterile swabs and sterile saline solution

Treatment cycles:

- A sustained treatment success can be expected only after repeated treatments.
- Usually 3-5 treatments are carried out, but depending on scar depth and structure also more treatments may be necessary.
- The interval between two treatments should be at least 4-6 weeks to make sure that the wound healing cascade is completed.

Contraindications:

- Treatment of mucous membranes and the eye - the skin around the eyes should be treated with extreme caution
- Patients with or a history of keloids
- Treatment of skin areas with skin cancer, active herpes simplex, active skin infections, raised moles, warts, solar keratoses

Safety notes

- Though no clear contraindication for microneedling, special attention is advised in patients with increased bleeding tendency or wound healing problems, e.g. due to uncontrolled diabetes mellitus or anti-coagulant therapy, as well as patients receiving chemotherapy, radiotherapy or high doses of corticosteroids
- Pregnant / nursing women: since there are no clinical results on pregnant / nursing women, we recommend to postpone the treatment
- Children / adolescents: the clinical study results were obtained mainly in adult volunteers. Thus, the treatment is generally recommended to be used in patients "over 18 years". However, there is no evidence that the treatment in adolescents and children constitutes an increased risk to the safety of the patient.
- When combining the microneedling treatment with other therapies, the other manufacturer's instructions for use have to be followed

Accessory symptoms of the treatment:

Transient redness and swelling of the treated skin, pain, microbleedings for minutes

Potential adverse effects:

Skin infections due to inadequate disinfection, inflammatory response (caused by pyrogens), transient pigmentation disorders, scarring when a tendency to keloid formation is present, granuloma after the use of active substances during the microneedling treatment.

A good skin hygiene before and after the procedure is advised to further minimise the risk of infection.

Prolonged sun or other UV exposure is not recommended after the procedure for at least 8 days to minimize the risk of hyperpigmentation. Patients should use a sun screen that provides a high level of UVA and UVB protection.

The stainless steel of the Dermaroller microneedles contains ca. 8% nickel. Due to the short treatment duration, it is very unlikely that a patient develops an allergic reaction.

Disposal

Used needle modules pose a risk of injury and a potential health hazard due to contact with blood. They are classified as „sharps“ (European Waste Catalogue 18 01 01) and must be collected in separate, closed containers and then disposed of safely.

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Révision 2018.11.06

 **DERMAROLLER GmbH**
Wilhelm-Mast-Str. 8
38304 Wolfenbuettel
Allemagne
Tel. +49 5331 7108480
info@dermaroller.de

MODULES D'AIGUILLES: NM6xx, NM6xx LS

Dispositif médical selon la directive 93/42/CEE – MADE IN GERMANY

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