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
CE
0482

Hersteller

Dermaroller GmbH
Wilhelm-Mast-Straße 8
38304 Wolfenbüttel
Germany

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
Manufacturer
Dermaroller GmbH
Wilhelm-Mast-Straße 8
38304 Wolfenbüttel
Germany


TABLE OF CONTENTS

1.	General notes	5
1.1	Operating manual.....	5
1.2	Representation of the warnings in this manual	5
1.3	Basic principles / Disclaimer	5
1.4	Changes to the device and Technical progress.....	6
2	Important safety instructions	6
2.1	Intended use	6
2.2	Indications.....	6
2.3	User groups	6
2.4	Patient group.....	7
2.5	Important hygiene and safety instructions	7
2.6	Contraindications and adverse effects	9
3	Delivery contents, accessories, spare parts, consumables	10
3.1	Delivery contents	10
3.2	Accessories, spare parts, consumables	11
4	Commissioning, operation and decommissioning	12
4.1	Control unit	12
4.2	eDermastamp Needle Module (consumable)	13
4.3	Handpiece	14
5	Device care, cleaning, and disinfection	15
5.1	Device care instructions	15
6	Device Maintenance	16
6.1	Safety inspection.....	16
7	Functional failure, repairs, complaints, disposal	16
7.1	Functional failure / repairs.....	16
7.2	Complaints.....	16
7.3	Disposal.....	17
8	Technical parameters	17
9	Operation, transport and storage conditions	17
10	Manufacturer Declarations	17
10.1	Warranty statement.....	17
10.2	Electromagnetic compatibility (EMC)	18
10.3	Manufacturer's declarations on electromagnetic emission	18
10.4	Manufacturer's declarations on electromagnetic immunity.....	19
10.5	Recommended safe distances from portable and mobile RF communication devices	21
10.6	List of symbols used on product and packaging	22
11	Notes	24

1 General notes

1.1 Operating manual

This operating manual explains the safe and proper use of the device and serves to meet health and safety requirements. All persons who implement, use, maintain or check the device must read and observe these instructions before use.

This operating manual is part of the device and must be at the operator's disposal at all times. If the device is passed on to a third party, this manual must also be passed on.

For the safe operation of the device, please also note

- the instructions for use of the eDermastamp Needle Module
- the Safety Data Sheets of disinfectants and cleaning agents
- the legal provisions on safety at work and microneedling

1.2 Representation of the warnings in this manual

	General Warning Sign – identifies risk of injury
	indicates a hazard that can result in death or serious injury if the hazard is not avoided
	indicates a hazard that may result in minor to moderate injury if the hazard is not avoided
	indicates potential risks that could result in damage to the environment, property or equipment if the hazard is not avoided
	Sharps hazard

1.3 Basic principles /Disclaimer

This device may only be operated by authorized, competent and qualified personnel and only for the intended uses stated in this operating manual (see chapter 2.2 "Indications")! Details can be found in chapter 2.3 "User groups".

These operating instructions do not replace the essential requirement for proper training in microneedling. A training course will provide information regarding any possible hazards and risks that are potentially associated with the treatment.

The user is obliged to adhere to the principles of cleanliness and hygiene and to work only with sterile or disinfected working materials. This explicitly includes the single use of sterile eDS needle modules. Further information can be found in chapter 2.5 "Important hygiene and safety instructions".

The contraindications in chapter 2.6 must be observed.

Inform your patients about any potential risk prior to the treatment.

The manufacturer accepts no liability for any adverse result, resulting from improper installation or use of the device, or for any damage resulting from this.

The use of the device is only intended in combination with the original eDermastamp Needle Modules from Dermaroller GmbH.

1.4 Changes to the device and Technical progress

This device may not be changed by the user in any way! Any change to the device involves the risk of serious injury or death, e.g. by an electrical shock.

The manufacturer reserves the right to alterations in terms of technical progress.

2 Important safety instructions

2.1 Intended use

The intended use for the system „eDermastamp and eDermastamp Needle Modules“ is

Minimal-invasive perforation of the epidermis / dermis (Microneedling)

The effect of a microneedling treatment with the eDermastamp system is puncturing the outer skin layers, epidermis and - depending on needle length - also the dermis, in a minimally invasive way without causing open wounds. The microchannels created by the needling will close up quickly again. The physiological repair process for renewal of injured skin tissue (the wound healing cascade) enables regeneration of epithelial function within a short time. Because this is connected to collagen production in the skin, this process is also called „percutaneous collagen induction“ (PCI) or "collagen induction therapy" (CIT).

2.2 Indications

Typical medical indications for Microneedling are

- *Treatment of disease and injury related skin lesions, especially atrophic and hypotrophic scars (such as acne scars), as well as mature hypertrophic scars on the whole body*
- *Treatment of Striae distensae*

In addition, the eDermastamp is suitable for perforating the skin surface to facilitate the introduction of active substances into the skin whose safety has been proven in suitable studies for the respective application.

2.3 User groups

The device may only be operated by the following persons:

- Healthcare professionals, such as physicians, qualified in applying the Microneedling technique

The qualification to carry out Microneedling treatments includes:

- Basic knowledge of the Microneedling procedure
- Knowledge of the hygiene and safety regulations
- Knowledge about contraindications, risks and possible side effects

2.4 Patient group

Treatment with the eDermastamp is suitable for patients of all skin types according to Fitzpatrick (I-VI) and all scar depths (I-IV). The available clinical study results were predominantly obtained with adult subjects. Therefore the age recommendation "over 18 years" is given for the treatment. However, there is no evidence that treatment in adolescents poses increased risks to patient safety. There are significant clinical data generated with the manually operated "Dermaroller MC" device by Dermaroller GmbH in the treatment of burn scars in children that demonstrate the safety of the treatment also in patients under the age of 18 years.

Parts of the body to which the eDermastamp is applied

The eDermastamp is generally applicable on almost all skin areas. Exceptions are described under 2.6 "Contraindications".

2.5 Important hygiene and safety instructions

A microneedling treatment can transmit infectious diseases if the required hygiene standards for invasive procedures are not met. Improper handling can result in injury to the user or patient. Therefore, the following hygiene and safety instructions must be strictly followed:

2.5.1 Hygiene instructions

- The microneedling procedure must be carried out in a hygienically clean room suitable for medical treatment under sufficient illumination
- Before and after each use, the device must be cleaned and disinfected according to the instructions in Chapter 5
- Contamination of the handpiece by back-flowing liquids is excluded by the patented design of the needle module. Should contamination nevertheless occur, operation must be interrupted immediately and adequate cleaning / disinfection carried out.
- The operator must wear a fresh pair of sterile gloves for each treatment to prevent cross-contamination.
- Before each procedure the patient skin has to be cleaned with a mild skin disinfectant. Choice of disinfectant must follow relevant national guidelines.
- Avoid touching freshly needled skin in general.
- Needle modules are sterile consumables for **single use only!**
- Needle Modules have to be used in sterile condition only! Sterility of needle modules is only guaranteed in an undamaged sterile packaging.
- When the Needle Module is attached to the handpiece the operator has to be careful that the Needle Module tip does not come into contact with surfaces that are insterile or have not been disinfected thoroughly, like clothes. If a

Needle Module was in contact other than with the disinfected skin area or a sterile glove it has to be replaced immediately.

- **ATTENTION:** a needle prick injury from a used or contaminated module might cause transmission of hazardous diseases. In case of a needle prick injury you should seek the advice of a medical doctor immediately.
- All consumables must be used up by the given expiry date. Adequate storage conditions have to be maintained.
- After the expiry date or after use, dispose of needle modules in appropriate special containers (Safety Box) as hospital waste. Local disposal regulations must be observed.
- Do not put any used modules in the household garbage (danger of infection)! Unused needle modules may be disposed of with household waste.

2.5.2 Precautions during operation

- The device system „eDermastamp + Needle Module“ may only be used by trained specialist personnel (details in section 2.3 "User groups").
- Check the device and accessories for damage before each use!
- In the event of (suspected) damage, never try to repair the device yourself! Only an authorized dealer may carry out / provide a repair.

Also, do not remove covers from the control unit, handpiece, power supply, or footswitch - there is a risk of electric shock.

- The eDermastamp may only be used with the associated authorized accessories (such as foot switch, adapters) and consumables (needle modules); Failure to oblige may harm the user or the patient!

- **NEVER** touch the needle module tip while the handpiece is in operating mode - sharps hazard!
- In non-operating state all needles must be fully retracted into the module tip. Defective modules must be exchanged carefully. Mind the injury hazard!

Restriction of the function due to electromagnetic interference

Mobile communication devices (mobile phones, WLAN routers, ...) can impair the correct functioning of the device by emitting electromagnetic radiation.

- Always set up the device in such a way that the protection distances to potential sources of interference recommended in section 10.5 are observed
- Do not use mobile / portable communication devices near the device
- Also inform patients of the possibility of such disorders

2.6 Contraindications and adverse effects

Subsequent contraindications and adverse effects are the result of careful analysis of clinical literature and clinical studies on microneedling. If the responsible operator has any doubt that the safety of the person to be treated cannot be guaranteed, e.g. due to accompanying diseases, the microneedling procedure has to be omitted or stopped immediately.

2.6.1 Contraindications

In the subsequent cases microneedling is contraindicated:

- treatment of mucosa or the eye
- treatment of skin materials with skin cancer
- patients with keloids or with a disposition for this

For the following patients, for reasons of caution, microneedling should be carried out only on the express advice of a treating physician:

- Patients with active herpes simplex, skin infections, raised liver marks, warts, photokeratosis (excluding the targeted treatment of photokeratoses in the context of photodynamic therapy)
- Patients with the disposition to bleeding or disturbed wound healing, i.e. caused by uncontrolled Diabetes mellitus
- Patients under anti-coagulant-therapy, chemotherapy, radiation therapy, high doses of cortisol.
- Patients planning for pregnancy, pregnant and lactating women
- Children and youths: Existing clinical study results were mostly achieved with adults. Therefore it is recommended to do treatments only on adult persons „over 18 years“. However there has been no indication that the treatment of youths includes a higher risk for patient safety.

The use of topical products / active ingredients in conjunction with Microneedling is only permitted if the safety of these products has been proven for the respective application in suitable studies.

2.6.2 Potential adverse effects

Compared with other methods of scar / skin treatment, e.g. by laser light, side effects caused by Microneedling are rare and usually mild.

Known transient reactions to microneedling applications with the Needle module

- Pain, redness, swelling of the treated skin areas
- Micro-bleeding („blood dew“) for minutes

Adverse reactions

- Infection after insufficient skin disinfection
- Inflammatory reaction (caused by pyrogens)
- Transient pigmentary abnormalities
- Scar formation with predisposition for keloids
- Granuloma formation after introduction of substances during microneedling

All newly treated skin should be protected from UV and sunlight at least 8 days, e.g. by using a sunblock with a light protective factor of 50.

3 Delivery contents, accessories, spare parts, consumables

3.1 Delivery contents

Control unit, power supply with 4 power supply adapters (EU, UK, AUS, USA), anodized base plate for control unit and handpiece, handpiece and connection cable, footswitch with connection cable, Operating instructions, safety inspection sheet

The product life envisaged by the manufacturer covers 10 years for all parts included in the delivery contents, provided that all operating, storage, and transport conditions as well as all measures for maintenance / safety inspection described in this manual have been observed.

Due to device technology advancements, a replacement of damaged parts can not be guaranteed after this time.

The consumables "eDermastamp Needle Modules" are not included in the delivery contents. They can be ordered through the authorized eDermastamp dealer.

3.2 Accessories, spare parts, consumables

Art.No.	Article	Units / PU
AD-1.21-MN-MP	e-Dermastamp	1
MAD-1.2-MN-MP	Anodized Base Plate	1
CME-2.0-MP	Handpiece eDS 1b	1
5E-G765	Pin cable for Handpiece	1
E-1170	Power supply GTM 41060-2515A 15 V	1
E-1153/-1154/-1155/-1157	Adapter USA/EU/UK/AUS	1 each
E-1017-ED-MP	Foot Switch	1
7E-AD1.21-MP -xx	Instructions for Use	1 per language
NM 615 LS	Needle module 1.5 mm	5

For all connections, only the supplied cables may be used.


4 Commissioning, operation and decommissioning

4.1 Control unit

Front view with operating and control elements

Rear view with connecting sockets

4.1.1 Commissioning

- Before first start-up, do a visual control to detect any transport damage.
- Use the device indoors only at a room temperature of 10 °C - 35 °C.
- Prior to using the device, please make sure that it has adapted to room temperature (a temperature differential of 10 °C requires at minimum three hours); otherwise, condensation could damage the electronics.
- Always place the device on a firm and level surface (table unit)
- The connection sockets for power supply unit and handpiece as well as the on / off switch are located at the back of the device and must always be accessible.
- The device may only be operated with the power supply stated on its nameplate. Line voltage must match the primary voltage imprinted on the power supply!
- Insert the country-specific adapter into the supplied power supply unit by pulling down the part of the power supply marked with an arrow and clicking on the corresponding adapter.
- The device is ready for operation when it is connected to the power supply. This is indicated by an LED next to the word "Power".

CAUTION Safe operation cannot be guaranteed:

- if the device or its components show any sign of damage,
- if the device does not work properly,
- after longer-term storage outside of the defined temperature range
- in case of transport damage if the device has been improperly transported.

In these cases the device should always be returned to the authorized specialist dealer/ supplier for inspection.

4.1.2 Operation

- Switch on the control unit by pressing the power switch button on the rear side of the device. The frequency indicator lights up on the display.
- Put the handpiece into operation by pressing the "GO" button.
- To switch off, press the "GO" button again.
- You can connect up to two handpieces to the device. Select the desired handpiece with the keys "1" or "2"
- You can also operate the unit with the supplied footswitch. The handpiece is switched on and off by stepping on the footswitch.

4.1.3 Setting the puncturing frequency

- Set the puncturing frequency using the arrow keys "up" (faster) / "down" (slower) when the device is switched on. The selected frequency is shown in the display.
- Microneedling is recommended at a puncturing frequency of at least 100 Hz to minimize the risk of "ripping" at the needle insertion point.
- Start with the preset frequency of 100 Hz and then adjust the frequency depending on skin condition and working speed. Too high stitching frequencies can lead to injuries to the patient.

4.1.4 Decommissioning


- Disconnect the unit by unplugging the power adapter cable from the unit.
- The device may only be transported or stored in its original or an equivalent packaging.

4.2 eDermastamp Needle Module (consumable)

- Before using the needle module, it is essential to observe chapter 2.5 "Important hygiene and safety instructions" and the Needle Modules' operating instructions!
- Only use the eDermastamp device with the *sterile original eDermastamp Needle Modules* by Dermaroller GmbH and for the indications given in chapter 2.2.
- The needle modules are single use items. They are delivered in sterile condition and individually packed in a sealed sterile pouch.
- The manufacturer guarantees sterility of the modules up to their attached expiry date, if kept in the original sealed and undamaged sterility packaging and under the specified transport and storage conditions (indicated on the modules' label).
- Store the needle modules in their original packaging until use. The Needle Module should be taken from the sterile pouch immediately before use.
- The label of the needle modules contains the lot number, date of manufacture and expiry date. The lot number is used to uniquely identify the production batch and must be noted down in the patient record. This number must be communicated to the manufacturer in case of complaint.

- Only needle modules of type NM 6xx LS (e.g. NM615 LS) may be used with the hand piece „CME-2.0-MP“.

4.2.1 Mounting the Needle Module to the handpiece

- Directly before use, take the Needle Module from the sterile packaging.
- Screw it clockwise onto the handpiece.
A noticeable resistance marks fixation of the module in the handpiece

4.2.2 Removing the Needle Module from the handpiece

Before removing the needle module, make sure that all needles are completely retracted into the module tip to prevent injury.

- To remove, unscrew the needle module counterclockwise.
- Make sure that no remaining liquid runs back into the handpiece by holding the handpiece appropriately.

4.2.3 Disposal of the Needle Module

- Never dispose of used needle modules in household waste (danger of infection through the needles), but in suitable special containers (Safety Box) as hospital waste. Observe the local disposal regulations.

4.3 Handpiece


The handpiece, in particular the drive located therein, is a mechanical precision instrument and must not be soiled. Therefore, no liquid may enter the drive.

Damage due to penetration of substances into the drive shaft shall not be recognized by the manufacturer as a reason for complaint.

Care must also be taken to ensure that the cable connecting the handpiece is not subjected to tensile stress or kinking.

4.3.1 Adjusting penetration depth

- Setting of the penetration depth must be carried *on the running handpiece*.
- Twist the upper part of the handpiece against the lower one (see illustration). Do not touch the tip of the needle module!

- The scale on the handpiece indicates the approximate needle length (0 to ~ 1.5 mm). The optimal needle protrusion differs depending on the type of application, the skin area to be treated and other factors.
- For treatment of acne scars, it is recommended to set the needle length to 1.5 mm.

5 Device care, cleaning, and disinfection

- Protect the handpiece from damage by placing it in on the handpiece holder on the base plate after each use.

5.1 Device care instructions

Always disconnect the device from the main circuit before starting any cleaning and maintenance work!

- Before and after each use of the device, clean and disinfect the surfaces of control unit, power supply, handpiece holder, and the handpiece with a mild disinfectant (such as 50% isopropyl alcohol) and a soft cloth.
- Use only approved cleansing agents / disinfectants which are compatible with the materials and surfaces of the device parts. Section 5.1.1 and 5.1.2 will help you to choose appropriate agents. Follow the applicable national regulations.

5.1.1 Recommended disinfectants

Supplier	Product	Exposure time
Antiseptica	Big Spray „new“	1 - 5 min
Bode Chemie	Bacillol	30 s - 1 min
Ecolab	Incidin Foam	1 - 2 min
Schülke & Mayr	Mikrozid Liquid	1 - 2 min

The hand piece must never be immersed in disinfectant as otherwise the parts inside might be damaged. Such damage is not covered by the manufacturer's warranty.

5.1.2 Material resistance

Resistant against weak acids (e.g. boric acid, acetic acid $\leq 10\%$, citric acid $\leq 10\%$), aliphatic hydrocarbons (e.g. pentane, hexane), ethanol as well as the most common inorganic salts and their aqueous solutions (e.g. sodium chloride, calcium chloride, magnesium sulfate).

Not resistant against strong acids (e.g. hydrochloric acid $\geq 20\%$, sulphuric acid $\geq 50\%$, nitric acid $\geq 15\%$), oxidising acids (e.g. peracetic acid), bases (e.g. sodium hydroxide solution, ammonia; all substances with a pH-value >7), aromatic/halogenated hydrocarbons (e.g. phenol, chloroform) as well as acetone und benzine.

6 Device Maintenance

6.1 Safety inspection

- The device has passed the technical safety inspection after manufacturing. The inspection protocol is included in the delivery contents.
- **Every 24 months** the device **must** be subjected to a renewed safety inspection. Hand over the device including hand piece, foot switch, power supply, and cables to your specialist dealer for distribution to an authorised servicing partner.
- Before each use, make sure that the date specified on the safety inspection sticker attached to the device is not exceeded

7 Functional failure, repairs, complaints, disposal

7.1 Functional failure / repairs

- If the device does not work properly, first disconnect all components such as power supply and handpiece from the device.
- Check all connections and reconnect all components; the device should now work properly again.
- If malfunction cannot be eliminated, hand over the device to your authorized dealer.
- Do not undertake any repair attempts yourself!
- If the device can not longer be operated in a safe manner, it is your duty to take it out of service and safeguard it against unsupervised or unauthorised use.

7.2 Complaints

- Complaints are only to be handled by your authorized dealer/supplier.

7.3 Disposal

According to the German „Electrical and Electronic Equipment Act“ (ElektroG), which puts EU-Directive 2012/19 (“WEEE Directive”) into German law, all manufacturers are obliged to take back electronic waste resulting from *commercial use* of the equipment and to dispose of it in a professional manner.

As manufacturer of the eDermastamp, we are registered with the German „Stiftung Elektro-Altgeräte-Register“ (WEEE registration number 18153250). If you need to dispose of your eDermastamp device or electronic accessory, please contact us. The device must not be disposed of with household waste or in public collection points for electronic waste.

8 Technical parameters

Type	AD-1.21-MN-MP
Rated voltage	15 V – (DC)
Power consumption	7 VA max.
Power supply	GTM41060-2515A
Protection class	2
Puncturing frequency	100 to 150 Hz ± 10%
Drive	Precision engine-DC
Operation mode	Continuous operation
Width x height x depth	230 x 45 x 185 mm
Weight of the hand piece incl. cable	ca. 100 g
Overall weight of the device	ca. 1100 g
Medical device classification	Class IIa

9 Operation, transport and storage conditions

Operation conditions	
Ambient temperature	+10 °C to +35 °C
Relative humidity	30% to 75%
Transport and storage conditions	
Protect from moisture and direct sunlight/heat	

10 Manufacturer Declarations

10.1 Warranty statement

With this device you have acquired a high quality medical device. The reliability of the device is ensured by the latest test methods. The respective national legal period of manufacturer warranty and liability is applicable for device defects which are caused by material or processing defects. No liability is assumed for consequential damages. For damages which are due to improper handling of the device or failure to follow our instructions, no warranty is accepted.

10.2 Electromagnetic compatibility (EMC)

Medical electrical equipment is subject to special precautions with regard to electromagnetic compatibility (EMC). This device may only be used for the purpose described in this manual and installed and commissioned in accordance with the EMC guidelines in this manual. Portable and mobile RF communication devices such as mobile phones can affect medical electrical equipment.

The DIN EN 60601-1-2: 2015 / IEC 60601-1-2: 2014 standard was implemented with the aim of ensuring that EMC requirements are met and insecure device statuses are avoided. The standard defines the level of electromagnetic disturbances for medical devices. The eDermastamp Microneedling devices are compliant to the standard in terms of both immunity and emissions.

The manufacturer guarantees conformity of the device with the EMC requirements only when original accessories and spare parts are used. Using other accessories may result in increased electromagnetic emissions or reduced immunity to electromagnetic interference.

10.3 Manufacturer's declarations on electromagnetic emission

The device is designed for use in an electromagnetic environment as specified below. The user of the device should ensure that it is operated in such an environment

Emission measurements	Conformity	Electromagnetic ambient conditions - guideline
RF emission as per CISPR 11	Group 1	The device uses RF energy only for its internal function. Therefore its RF emissions are very low and an influence on nearby electronic devices is unlikely.
RF emissions as per CISPR 11	Class B	The device is designed for use in all facilities, including residential areas and areas directly connected to a public power supply network that also supplies buildings used for residential purposes.
Harmonics as per IEC 61000-3-2	Not applicable	
Voltage fluctuations / flicker as per IEC 61000-3-3	Not applicable	

10.4 Manufacturer's declarations on electromagnetic immunity

The device is designed for use in an electromagnetic environment as specified below. The user of the device should ensure that it is operated in such an environment

Immunity tests	IEC 60601 test level	Compliance level	Electromagnetic environment - guidelines
Electrostatic discharge (ESD) as per IEC 61000-4-2	± 6 kV Contact discharge ± 8 kV Contact discharge	Compliant	Floors should be covered with ceramic tiles. If the floor is covered with a synthetic material, the relative air humidity should be at least 30%.
Electrical fast transient disturbances /bursts as per IEC 61000-4-4	± 2kV for power lines ± 1 kV for input / output lines	Compliant	The quality of the supply voltage should be equivalent to that for a typical business or hospital environment.
Surges as per IEC61000-4-5	± 1 kV Differential mode voltage ± 2 kV Common mode voltage	Compliant	The quality of the supply voltage should be equivalent to that for a typical business or hospital environment.
Voltage dips, short interruptions and fluctuations in supply voltage as per IEC61000-4-11	< 5% UT (> 95% dip of UT) for ½ period 40% UT (60% dip of UT) for 5 periods 70% UT (30% dip of UT) for 25 periods < 5% UT (> 95% dip of UT) for 5 s	Compliant	The quality of the supply voltage should be equivalent to that for a typical business or hospital environment. If the user of the device requires continued operation even in the event of interruptions to the power supply, it is recommended to supply the device from an uninterruptible power supply or battery.
Magnetic field at supply frequency (50/60 Hz) as per IEC 61000-4-8:	3 A/m	Compliant	Magnetic fields at the mains frequency should be equivalent to the typical values found in the business and hospital environment.

* UT is the alternating mains voltage prior to application of test levels

Immunity tests	IEC 60601 test level	Compliance level	Electromagnetic environment - guidelines
			Portable and mobile radio equipment is not used at a distance from the device, including its cables, that is less than the recommended safe distance, which is calculated using the appropriate equation for the transmission frequency. Recommended safe distance:
Conducted RF disturbances as per IEC 61000-4-6	3 V _{eff} 150 kHz to 80 MHz	Compliant	$d = 1,2 \sqrt{P}^{**}$
Radiated RF disturbances as per IEC 61000-4-3	3 V/m 80 MHz to 2,5GHz		$d = 1,2 \sqrt{P}^{**}$ 80 MHz to 800 MHz $d = 2,3 \sqrt{P}^{**}$ 800 MHz to 2,5 GHz
<p>** where P is the nominal power of the transmitter in watts (W) as specified by the transmitter manufacturer and d is the recommended safe distance in meters (m). At 80 MHz and 800 MHz, the higher value applies.</p> <p>These guidelines may not apply in all situations. The propagation of electromagnetic waves is influenced by absorption and reflection from buildings, objects and people.</p>			
			<p>The field strength of stationary radio transmitters at all frequencies in accordance with a test on site ^{a)} should be less than the compliance level ^{b)}.</p> <p>Interference is possible in the vicinity of devices marked with the symbol shown on the left.</p>
<p>^{a)} It is not possible to predict theoretically the field strength of stationary transmitters (such as base stations for mobile phones or land mobile services, amateur radio, AM and FM radio) with accuracy. To determine the electromagnetic environment with respect to stationary transmitters, a study of the location should be considered. If the measured field strength at the device location exceeds the above mentioned compliance level, the device must be observed during normal operation to prove that it is working as intended. If any unusual performance characteristics are observed, additional measures may be necessary as, for instance, a change of orientation or relocation of the device.</p> <p>^{b)} Above the frequency range of 150 kHz to 80 MHz, the field strength is less than 3 V/m.</p>			

10.5 Recommended safe distances from portable and mobile RF communication devices

The device is designed to be used in an electromagnetic environment in which RF disturbances are controlled. The user of the device can prevent electromagnetic interference by keeping a minimum distance between the device and portable and mobile RF communication devices (transmitters) as indicated below, depending on the power of the communication device.


Nominal power of transmitter [W]	Safe distance depending on transmitting frequency [m]		
	150 kHz to 80 MHz $d = 1,2 \sqrt{P}$	80 MHz to 800 MHz $d = 1,2 \sqrt{P}$	800 MHz to 2,5 GHz $d = 2,3 \sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

If the maximum output power for a transmitter is not given in the above table, the recommended safe distance d in metres can be determined by using the equation given in the respective column, where P is the maximum output power of the transmitter in Watt (W), as specified by the transmitter manufacturer.

At 80 MHz and 800 MHz the higher frequency range is applied.

These guidelines may not apply in all situations. The propagation of electromagnetic waves is influenced by absorption and reflection of buildings, objects, and people.

10.6 List of symbols used on product and packaging



IP21	Housing protects against coarse dust and water droplets.
	Applied part type B; Application part provides protection against electric shock and leak current.
	DC connection / inner pin positive
	The device complies with the requirements of the Directive 93/42/EEC concerning medical devices.
	„Refer to instruction manual“
	„Manufacturer“
	„Manufacturing date“
	„Catalogue number“
	„Serial number“
	„Lot (batch) code“
	„Sterilized with ethylene oxide“
	Expiry date
	„Attention! Read warning notes!“ / General Warning
	„Keep in a dry place“
	„Fragile - handle with care “
	„Product shall not be disposed of in household waste!“
	„Do not use if packaging is damaged!“
	„Do not reuse!“
	„Do not resterilise!“
	Control sticker from safety inspection

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