

# Improved Handpiece Connection

The handpiece and the control unit are now permanently connected. The previous plug connection has been removed: the device is easier to clean, and the plug can no longer be fitted incorrectly by the user. For servicing, the handpiece can still be replaced as before.

## BEFORE



**Previous version.** Handpiece connected to the control unit via a plug.

## NEW



**New version.** One continuous cable, permanently attached to the handpiece.

## Development objectives

- 01** Permanent connection between handpiece and control unit
- 02** Easier cleaning
- 03** Fewer interfaces
- 04** Simple replacement by the service technician in case of repair

## Benefits of the new version

- ✓ **Smooth cable** that is easy to clean
- ✓ **No more gaps** at the plug that are hard to clean
- ✓ **Handpiece replacement** by the service technician still possible
- ✓ **Incorrect handling by the user is ruled out**
- ✓ **Cleaner, symmetrical design**
- ✓ **Existing packaging** can still be used



**BEFORE** Connection at the rear of the device.



**NEW** Connection at the front.



**UNCHANGED** Packaging and accessories.